




# FLOWGUIDE

INCREASED PERFUSION DURING DEEP INTUBATION

# FLOWGUIDE

INCREASED PERFUSION DURING DEEP INTUBATION

ENABLES CONVENIENT USE WITH OTW DEVICES



DESIGNED FOR THE  
TRANSRADIAL  
APPROACH



| MODEL | REQUIRED GUIDE CATHETER ID     | FLOWGUIDE ID     | RAPID EXCHANGE LENGTH | USABLE LENGTH |
|-------|--------------------------------|------------------|-----------------------|---------------|
| FG60F | 6F ID $\geq 0.070$ " (1.78 mm) | 0.056" (1.42 mm) | 15 cm                 | 150 cm        |

Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.  
Caution: Federal Law (U.S.) restricts this device to sale by or on the order of a physician.  
FlowGuide™ is a trademark of IMDS.

© IMDS 2024 ITM220037 R3.0

[WWW.IMDS.NL](http://WWW.IMDS.NL)

Ceintuurbaan Noord 150  
9301 NZ Roden  
The Netherlands  
[sales@imds.nl](mailto:sales@imds.nl)

Phone +31(0)50 8200230  
Fax +31(0)50 8200231

