


Finethrough[®]

Cardiovascular Guidewire

SINA



Parameters

O.D.: 0.014"/Max0.36mm

Total length: 180cm/300cm

Tip radiopacity: 3.0cm

Tip Shape: Straight / J

***FEATURES**

SS Core in one piece with PTFE coating

Providing improved torque transmission and guide wire trackability.

Superior Hydrophilic coating


Enables excellent controllability during complex lesions.

Distal dual coil design

Ensures beyond expectations of tip durability and torque performance.

Proximal marker


For clear identification of Tip Load



- Tip load: 0.7gf, 0.5gf
- Tip radiopacity: 3cm
- Spring coil length: 25cm
- Hydrophilic coating length: 28cm
- Proximal coating: PTFE

SS304V Core to tip design

One-piece core wire supports the entire guide wire from the proximal to the distal end. This design transmits the guide wire torque from one end to the other without loss.



ORDERING INFORMATION

REF.NO.	Diameter	Length	Tip Load	Spring coil length	Radiopaque length	Tip Shape
NS1403120AS	0.014 inch	180cm	0.7gf	25cm	3cm	Straight
NS1403120AJ	0.014 inch	180cm	0.7gf	25cm	3cm	J
NS1403120CS	0.014 inch	300cm	0.7gf	25cm	3cm	Straight
NS1403120CJ	0.014 inch	300cm	0.7gf	25cm	3cm	J
NS1403110AS	0.014 inch	180cm	0.5gf	25cm	3cm	Straight
NS1403110AJ	0.014 inch	180cm	0.5gf	25cm	3cm	J
NS1403110CS	0.014 inch	300cm	0.5gf	25cm	3cm	Straight
NS1403110CJ	0.014 inch	300cm	0.5gf	25cm	3cm	J



Add

No.925 Jinyuan yi Road,
Shanghai, P.R.C.



TEL

021-59140046



FAX

021-69000013



P.C.

201803



WEB

www.kdl-int.com



Shanghai Kindly Medical Instruments Co., Ltd