


RECROSS



SCAN CODE FOR
INSTRUCTION VIDEOS

MODEL	GUIDE CATHETER COMPATIBILITY	DISTAL SHAFT PROFILE	PROXIMAL SHAFT PROFILE	GUIDEWIRE COMPATIBILITY	USABLE LENGTH
RCM4025	SF	2.3F/3.3F	2.6F/3.4F	0.034"	140 cm

Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.
Caution: Federal Law (U.S.) restricts this device to sale by or on the order of a physician.
In-vitro test data on file at IMDS.
Bench test data may not be indicative of clinical performance.
ReCross™ is a trademark of IMDS.


© IMDS 2023

ITM200184 R2.0




Centaurbaan Noord 150
9501 NZ Breda
The Netherlands
Sales@imds.nl
+31(0)40 8200050
+31(0)50 8200221

TOOLS FOR EVERYDAY EFFICIENCY IN COMPLETE REVASCULARISATION




WWW.IMDS.NL

PRECISE TARGETING OF SIDE BRANCH
OCCLUSION



STRONG FRONTLINE CROSSING AND
PENETRATION WHEN NEEDED



ANTEGRADE
DUAL ACCESS (ADA)

RECCROSS

UNPARALLELED GUIDEWIRE REDIRECTION OPTIONS

360° CIRCUMFERENTIAL ACCESS DUAL LUMEN
ENABLES UNIQUE & EFFECTIVE PROCEDURES

OVAL DISTAL SHAFT
MINIMIZING CROSSING PROFILE



RECCROSS


VERSATILITY FOR **EVERYDAY** EFFICIENCY

ENHANCER Rx



6F GUIDING CATHETER
COMPATIBILITY

SUBINTIMAL GUIDEWIRE
REDIRECTION (SGR)



REMOVABLE STYLET
PROVIDES SIGNIFICANT PUSHABILITY IMPROVEMENT

HMATOMA
DECOMPRESSION