


SILENE

Peripheral Covered Stent™

Mono Layer Covered Stent



- A versatile stent to treat a wide range of vascular pathologies
- Clinically proven results in patients worldwide*
- "Even in the most difficult anatomical conditions we could appreciate the properties of the Silene stent. In particular, the low profile which allowed us easy access from the humeral site, the excellent navigability and the precision of release, the possibility of overdilation and high radial force."**



*Choose Reliability.
Choose Silene.*

Technical Data

Specifications	Description
Stent Design	Direct-Stent® - Patented Hybrid Cell Design
Stent Material	L605 Cobalt Chromium Alloy
Stent Covering Polymer	ePTFE Microporous Polymer, IND 90 - 120 µm
Wall Thickness	0.08 mm/0.003"
Balloon System	Over the Wire (OTW)
Balloon Material	Semi-Compliant – Nylon Blend
Optimal Deployment Pressure	9 ATM
**Nominal Pressure	7 ATM
Rated Burst Pressure	10 ATM
Shaft Size	5/6F
*Minimum Sheath Size	Ø6 - 7 mm: 6Fr; Ø8 - 9 mm: 8Fr; Ø10 - 12 mm: 9Fr
Pre-mounted Profile	1.7 - 3.0 mm
Tip Entry Profile	0.039"
Recommended Guidewire	0.035"
Delivery Catheter Length	~ 130 cm
Stent Length	14 - 70 mm
Stent Ends Flaring	6.7 >>> 9 mm; 8, 9, 10 >>> 14 mm; 12 >>> 16 mm
Radial Force	10 - 9 mm 2.0N Force

* See individual manufacturer for Sheath Size Fr equivalent

** Assure full deployment of stent. Deployment pressures should be based on lesion characteristics.

Ordering Information

Stent Diameter (mm)	Stent Length (mm)	Ordering Number
6.0	15, 25, 35, 50, 70	75 Stent Diameter (with no decimal point) Stent Length For example: Stent Diameter: 6.0 Stent Length: 15 -> Ordering number: 756015
7.0	15, 25, 35, 50, 70	
8.0	15, 24, 32, 50, 70	
9.0	15, 24, 32, 50, 70	
10.0	15, 23, 31, 48, 65	
12.0	14, 23, 31, 48, 65	

The content in this brochure is provided for informational purposes only and is not intended or recommended as a substitute for educational professional medical advice or diagnosis purposes. The information is for use only in countries with applicable health authority product registration.

INSITU TECHNOLOGIES® INC.
 539 Phalen Blvd, St. Paul, MN 55130, USA
 Tel: +1 651 389 1017
 Email: hello@insitu-tech.com
 Web: www.insitu-tech.com

FDA Registered Facility
 EN ISO 13485:2016 certified
 For International Distribution, Rev.09/01/2021
 ©2021 COPYRIGHT InSitu Technologies® Inc.



Silene™ Peripheral Covered Stent

Catalog Number	Description	Catalog Number	Description
756015	Silene™ 6.0 x 15	759015	Silene™ 9.0 x 15
756025	Silene™ 6.0 x 25	759024	Silene™ 9.0 x 24
756035	Silene™ 6.0 x 35	759032	Silene™ 9.0 x 32
756050	Silene™ 6.0 x 50	759050	Silene™ 9.0 x 50
756070	Silene™ 6.0 x 70	759070	Silene™ 9.0 x 70
757015	Silene™ 7.0 x 15	7510015	Silene™ 10.0 x 15
757025	Silene™ 7.0 x 25	7510023	Silene™ 10.0 x 23
757035	Silene™ 7.0 x 35	7510031	Silene™ 10.0 x 31
757050	Silene™ 7.0 x 50	7510048	Silene™ 10.0 x 48
757070	Silene™ 7.0 x 70	7510065	Silene™ 10.0 x 65
758015	Silene™ 8.0 x 15	7512014	Silene™ 12.0 x 14
758024	Silene™ 8.0 x 24	7512023	Silene™ 12.0 x 23
758032	Silene™ 8.0 x 32	7512031	Silene™ 12.0 x 31
758050	Silene™ 8.0 x 50	7512048	Silene™ 12.0 x 48
758070	Silene™ 8.0 x 70	7512065	Silene™ 12.0 x 65